

Assisi™ Animal Health



Portable-150™

Assisi Portable™-150

\$250

Assisi Portable™-150 is a non-invasive disposable targeted Pulsed Electromagnetic Field (tPEMF™) device weighing about two ounces. It has a 7.5 inch applicator. Treatment begins by simply placing the applicator over the treatment site and pressing the start button. The device will turn itself off once the 15-minute treatment is complete.

Features

- FDA Cleared
- Non-invasive and non-pharmacological
- Lightweight, conformable applicator is easy to position for therapy
- Treats directly through wraps, dressings and casts
- Easy to administer
- Painless therapy for patients
- 150 fifteen minute treatments
- To get 150 treatments from your Assisi Portable, please wait at least 2 hours between treatments. Treating more frequently will reduce the number of treatments from the device.



Portable Automatic™

Assisi Portable™ Automatic

Assisi Portable™ Automatic is a non-invasive disposable targeted Pulsed Electromagnetic Field (tPEMF™) device weighing about two ounces. It has a 7.5 inch applicator. Treatment begins by simply placing the applicator over the treatment site and pressing the start button. The Assisi Portable™ Automatic offers preprogrammed 15-minute treatments every 2 hours and can be integrated with post-surgical or wound care dressings. The device can be used manually or in the pre-programmed mode.

Features

- FDA Cleared
- Non-invasive and non-pharmacological
- Lightweight, conformable applicator is easy to position for therapy
- Treats directly through wraps, dressings and casts
- Easy to administer
- Treats manually or automatically
- Painless therapy for patients
- 100 fifteen minute treatments
- To get 100 treatments from your Assisi Portable, please wait at least 2 hours between treatments. Treating more frequently will reduce the number of treatments from the device.